



Three-in-One Bariatric Commode

ProBasics Bariatric 3-in-1 Commode has a wider seat and steel frame to support a higher weight capacity and is ideal for larger patients. It has three versatile functions; it can be used as a toilet safety frame, a commode, and a raised toilet seat. It comes complete with commode bucket and lid. The versatile design includes a rectangular pan and pan holder, armrests for added comfort, angled legs for stability and non-marring rubber tips.

FEATURES

- Can function as a free-standing commode, a toilet safety frame (with bucket removed) or as a raised toilet seat with arms.
- Durable steel-welded frame is attractive, strong and easy-to-clean
- Tool-free assembly, includes plastic snap on seat with lid
- The rectangular pan is accessible from both the front and back of commode
- Plastic arms provide added comfort
- Seat height readily adjusts in one-inch increments
- Color: gray frame, white accessories
- HCPCS Code: E0168

SPECIFICATIONS

Width b/w Arms:	21"
Outside Width:	30"
Seat Dimensions:	15.5"W x 15.5"D
Seat Opening:	8.25"W x 10"D
Adjustable Seat Height:	17" - 23"
Color:	Gray / White
Weight Capacity:	450 lb
Carton Dimensions:	29.5" x 17.5" x 15.5"
Shipping Weight:	35.5 lb
Warranty:	Limited lifetime on frame One year on parts

Item #	Description	UPC	UOM
BSB31C	ProBasics Bariatric 3-in-1 Commode	815067071852	2/cs



HEALTHCARE
HOME MEDICAL SUPPLY



Call or Visit Our Website



www.healthcaredme.com

For Further Details Please Contact Us.
Address: 2911 Carpenter Rd, Ann Arbor, MI 48108.
Phone: (734) 975-8668, Fax: (734) 975-6678
Email: info@healthcaredme.com

DmeHealthcare HealthcareHomeMedical