



## Bariatric Rollator with 8-inch Wheels

The ProBasics Bariatric Aluminum Rollator is made with an extra-strength, reinforced aluminum frame and has large reinforced forks on 8-inch wheels for improved maneuverability and stability. The padded seat and removable padded backrest provide a comfortable resting area. It features easy-to-operate ergonomic hand brakes, a water-resistant nylon pouch, and a built-in fold strap for easy storage and transport. Supports patient weights up to 400 pounds



Extra Wide Seat and  
Water-Resistant Pouch



Easy-to-Operate  
Ergonomic Hand Brakes



Reinforced forks and  
8-inch Wheels

### FEATURES

- Lightweight aluminum frame
- Two-button removable curved padded backrest
- Comfortable, extra wide flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 8" black casters with reinforced forks and soft grip tires, ideal for both indoor and outdoor use
- HCPCS\* Code: E0149

E0156 (seat)

### COLORS

Burgundy  
Finish



Blue  
Finish



### SPECIFICATIONS

Overall Width:	29"
Width b/w Hand Grips:	23.5"
Seat Dimensions:	19" (W) x 13" (D)
Seat-to-Floor Height:	22"
Weight Capacity:	400 lb
Adjustable Handle Height:	33.5" – 38.5"
Wheel Size:	8"
Product Weight:	18.7 lb
Carton Dimensions:	29.1" x 8.6" x 36.2"
Shipping Weight:	22.7 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLAB8BG	ProBasics Bariatric Aluminum Rollator	Burgundy	815067071746	1/cs
RLAB8BL	ProBasics Bariatric Aluminum Rollator	Blue	815067071739	1/cs



**HEALTHCARE**  
HOME MEDICAL SUPPLY



Call or Visit Our Website



www.healthcaredme.com

For Further Details Please Contact Us.  
Address: 2911 Carpenter Rd, Ann Arbor, MI 48108.  
Phone: (734) 975-6668, Fax: (734) 975-6678  
Email: info@healthcaredme.com

DmeHealthcare HealthcareHomeMedical