



Bariatric Transfer Bench

The ProBasics Bariatric Transfer Bench helps make transfers in and out of the bathtub safer and more comfortable. Constructed with an aluminum anodized frame, it features added cross braces for stability and a heavy duty blow-molded seat. It converts easily for left or right hand transfers. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Durable blow-molded bench and back have a textured surface and are contoured with drainage holes to reduce slipping
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Added cross brace for stability
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0248*

SPECIFICATIONS

Overall Footprint:	32.5"W x 21.5"D
Seat Dimensions:	25"W x 15.5"D
Adjustable Seat Height:	18" - 24.25"
Weight Capacity:	500 lb
Carton Dimensions:	33.5" x 20.5" x 17"
Shipping Weight:	28 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBTB	ProBasics Bariatric Transfer Bench	815067071975	2/cs