

JOBST® Elvarex® Soft Seamless



Now also
available in
Black!

The world's only range
of seamless flat-knit
compression garments.



#1 physician recommended



Comfort, Health and Style!

Optimum care with **seamless compression**


For over 50 years, JOBST® has been continuously using new knit technologies to improve the medical effectiveness of lymphatic compression therapy and to make compression garments more comfortable to wear.



This commitment to research and development has ensured that compression stockings made by BSN medical have become a key element in the treatment of lymphedema.

BSN medical introduced its JOBST® Elvarex® Soft Seamless range for the optimal compression treatment of lymphedema on the fingers and toes in 2007. These medical aids maintain their status as the world's first and only, tailor made, seamless flat-knit compression garments. There are no uncomfortable seams in the narrow spaces between fingers and toes, ensuring even and consistent compression.

The success story of JOBST® Elvarex® Soft Seamless: Research and development for the benefit of the patient.



JOBST® Elvarex® Soft Seamless

For hands and feet

Features and benefits:

JOBST® Elvarex® Soft Seamless is the brand which offers reliable help and relief in the case of lymphatic diseases of the hands and feet. Entirely seamless, this unique flat-knit range ensures outstanding comfort combined with the highest level of medical efficacy.

- **Made from a single piece of totally seamless material**

for uninterrupted compression, optimal edema control and outstanding comfort.

- **Anatomically shaped finger/toe separation**

makes the garment easy to put on and take off and provides an excellent fit, even in the critical area between the fingers and toes.

- **Flat edge finish at fingers and toes**

it is not bulky and prevents the ends from curling.

- **Special material and knit structure**

ensures a therapeutic micromassage effect.

- **Comfort zones at the edge**

easy to put on and comfortable to wear, especially in combination with compression sleeves and stockings.

- **Soft, durable and breathable knitted fabric**

helps improve patient compliance and enhances sensory perception.



JOBST® Elvarex® Soft Seamless

The seamless flat-knit compression glove

An overview of the technological innovations from BSN medical:



Now also available in Black!

The only flat-knit compression glove. Completely seamless. Extremely durable.

Special knitting technique in the palm of the hand creates anatomically correct shaping.


Double-covered inlay yarn in every stitch course ensures uninterrupted and controlled compression over the entire glove, even at critical areas.

Visually pleasing, breathable fabric is durable, yet soft to the touch and enhances sensory perception.

Seamless reverse edge. Constructed entirely without seams. This makes the glove particularly durable and wear-resistant.

Low profile, flat edge finish prevents curling.





Low profile, flat, seamless finger edges prevent curling.

Anatomically shaped finger gradient using 3D technique:
■ improves donning and removing

■ facilitates a particularly good fit and edema control in the spaces between the fingers

Perfect fit, even to the outer side of the hand, due to anatomically correct custom-made manufacture.

A further advantage of the seamless, flat-knit compression glove: The material and knitting structure of the glove give the advantage of a micro-massage effect exerted over the hand.

Comfort zone around the wrist makes the compression glove easy to put on and wear when worn in conjunction with the JOBST® Elvarex® armsleeve.

JOBST® Elvarex® Soft Seamless

The seamless flat-knit compression foot cap.

An overview of the technological innovations from BSN medical:



Comfort zone at the edge to promote wearing comfort when worn in conjunction with: Elvarex® and Elvarex® Soft leg compression garments.



Perfect fit, even to the outer side of the foot, due to anatomically correct custom-made manufacture.



Anatomically shaped toe gradient using 3D technique facilitates a particularly good fit and edema control in the spaces between the toes.



Low profile, seamless, flat elastic edges around the toes prevent curling.



Anatomically correct shaping due to a special knitting technique in the sole of the foot as well as the ball of the foot. Individual right and left foot cap for optimum fit.

Now also available in Black!

Pleasantly soft, visually pleasing, durable and breathable fabric. This low profile, flat-knit material is soft to the touch.

Double-covered inlay yarn ensures uninterrupted and controlled compression over the entire foot cap, even at critical areas.

Flat edge, low profile finish prevents curling.


Reverse edge constructed entirely without seams. This makes the foot cap particularly durable, wear-resistant and comfortable to wear.

A further advantage of the seamless, flat-knit compression foot cap: The material and knitting structure of the foot cap give the advantage of a micro-massage effect exerted over the foot.

Anatomically shaped opening for the small toe provides optimum comfort in wear.


JOBST® Elvarex® Soft Seamless

Seamless flat-knit compression garments for hands and feet.



- Perfect anatomical fit due to custom manufacture
- Now available in Beige and Black
- Individual right/left compression treatment for optimum fit
- Therapeutic efficacy confirmed by fit trials
- Manufactured from high-quality Nylon and Spandex fibers; latex-free

JOBST® Elvarex® Soft Seamless: the range



JOBST® Elvarex® Soft Seamless is only available as a custom-fit solution. The glove and the foot cap are available in standard beige and black color in compression classes 1 and 2. The AC¹ and AE solutions are available in the following versions: with fingers, with thumb and without fingers, with thumb. The foot cap includes an anatomically shaped opening for the small toe

For further information about JOBST® Elvarex® Soft Seamless, please do not hesitate to contact our Customer Service Department at (800) 537-1063.

Please refer to the product label and / or package insert for full instructions on the safe use of these products.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.jobst.com
To order toll-free: JOBST (+1) 800 537 1063

61340 RN

© 2012 BSN medical Inc.

1/1 physician recommended



Comfort, Health and Style!

REV 12/12