

State

Purchase Order No
Ordered by
Street
City

Fitter

PLEASE DIRECT ALL ORDERS TO

2920 Centennial Rd., Toledo OH 43617-1833 $\,$ toll-free:

MEDICAL, INC. www.gottfriedmedical.com www.gottfriedmedical.com <a href="mailto:volume="mailto


office: +1 419 474-2973 toll-free: +1 800 537-1968

Date	Fitter's Phone No.			
	Patient		Age	
	Ship to			
Zip	Street			
	Citv	State	Zip	

Α	Medial to Lateral (Inside to Outside)	E	
В	Anterior to Posterior (Front to Back)	Waist	
С	Inseam / Length of Stump	Last 1" fr	
E	Girth at Waist	C A Stump) End

Measurements (below) from Straight Edge of Blue Spine in: (check one) ☐ Inches ☐ Centimeters

Form F-032 - For Use With Paper Tapes



Please refer to our current price list for a description of these items

□ HOTLINE SERVICE

(Guaranteed delivery in 4-5 business days, or less) Additional Cost

Sex: ☐ M ☐ F

REQUIRED Allergies:

Allergic to silicone? ☐ Yes, ☐ No, ☐ NA | Allergic to metal? ☐ Yes, ☐ No, ☐ NA

WAIST LENGTH SUPPORTS - NOT SHOWN

See Form F-029

OPTIONS & EXTRAS:

- Zipper: ☐ Inside (or) ☐ Outside ☐ Open (or) ☐ Closed
- Garment Color No:
- Reinforcement: □ Knee Front | □ Knee Back
 - Inner Thigh
- Flex/Contracture Seam: ☐ Knee Back
- Elastic Band: Regular □ 1" □ 2" (or) Silicone □ 1" □ 2"
- Other Option(s):

INDICATIONS:

- 20-30 mm Hg: Varicose veins (mild), arterial insufficiency with venous insufficiency. Burns prevention of hypertrophic scars.
- □ 30-40 mm Hg: Varicose veins (moderate), assist fluid return, leg fatique, stasis dermatitis, postphlebitic syndrome, post surgical stripping of sclerosing, postfracture edema, prophylactic treatment of edema and phlebitis, lymphedema (moderate).
- □ 40-50 mm Hg: Chronic venous insufficiency, stasis dermatitis (severe), lymphedema (severe), chronic venous insufficiency (severe), orthostatic hypotension (moderate).

Contraindications: Non-ambulatory use, severe arterial insufficiency, cutaneous infection, acute hypodermatitis, wet dermatitis.