

14mm 16mm 18mm 20mm 22mm 24mm 26mm 28mm 30mm



**HEALTHCARE**  
HOME MEDICAL SUPPLY




# Flange Sizing Guide

Print this sheet at 100% custom scale for accuracy. Cut along the blue dotted lines to measure your nipple diameter.

## Where to Measure

Your flange size should be 3-5mm larger than your nipple diameter.



\*Be sure to measure BEFORE pumping.

## Helpful Hint

If you do not have access to a printer, use a coin to determine your nipple diameter.



32mm

34mm

36mm

If your nipple size is...

your flange should measure...

15-17mm

20mm

19-21mm


24mm

23-25mm

28mm

27-29mm

32mm



Call or Visit Our Website



[www.healthcaredme.com](http://www.healthcaredme.com)

For Further Details Please Contact Us.  
Address: 2911 Carpenter Rd, Ann Arbor, MI 48108.  
Phone: (734) 975-6668, Fax: (734) 975-6678  
Email: [info@healthcaredme.com](mailto:info@healthcaredme.com)

DmeHealthcare HealthcareHomeMedical